HYGROLABC1

SHORT INSTRUCTION MANUAL

Bench-Top Water Activity Indicator

General description

Congratulations on your purchase of the HYGROLAB C1 water activity indicator. Please read these short instructions carefully before installing the device. The detailed instruction manual can be found on the internet at: www.rotronic.com

The HygroLab C1 is a multifunction bench-top indicator that can read simultaneously up to 4 HygroClip 2 digital humidity-temperature probes. The HygroLab C1 is primarily designed for measuring the water activity (Aw) of foods, pharmaceuticals and cosmetics. The HygroLab C1 can also be used in other laboratory applications that require measuring humidity and temperature.


Power supply

The HygroLab C1 requires the power adapter model AC1211 (100...240 VAC – 12 VDC, 200 mA). As an alternative, the HygroLab C1 can also be powered via the USB port. In this case, the Ethernet port is inactive and cannot be used.

Probe inputs

The HygroLab C1 has 4 probe inputs designed for use with either HC2-AW or all standard HygroClip 2 probes.

Dimensions / Connections

HygroLab C1 back panel showing the four probe inputs as well as the USB and Ethernet ports:

Keypad

UP / DOWN

- Changes the probe being displayed
- When the menu is active, use to navigate the menu, make a selection or change a value

ON/OFF

Turns the instrument on or off

ENTER

- When the menu is active, use to confirm the selection of a menu item, effect a change of settings and confirm any change
- In the Aw mode, use to start, hold or stop the water activity measurement function

MENU


Activates the internal menu. Press this key again to go back to the previous sub-menu or to exit the menu

The ON/OFF key is used as follows:

- To turn the instrument on: press the ON/OFF key
- To turn the instrument off: press the ON/OFF key


Display

The LC display has a backlight which can be set to be on all the time or whenever a key is pressed. The backlight can also be disabled. Using the HygroLab C1 MENU > Display Settings > Mode, the display mode can be changed as shown below:

Standard:

- %RH
- Temperature
- Date and time

3-line display:

- %RH
- Temperature
- Calculated parameter
- No date and time

The display can also be configured to show a trend indicator on each line:

- ▲ increasing value
- ▼ decreasing value

- ◆ stable value

In the event of an alarm the symbol [!] appears to the right of the value

Frequently used settings

Unit system

Press the MENU key and select Device Settings > Unit Sys. Press ENTER to activate the Unit Sys menu item, use the UP or DOWN arrow key to change the unit system. Press ENTER to confirm and press MENU to exit.

The HW4 software can also be used to change the unit system.

WARNING: changing the unit system does not change the numerical value of the fixed barometric pressure used in some psychrometric calculations. Always verify and modify as required the numerical pressure value after changing the unit system.

Date and time

Press the MENU key and select Time Settings > Date or Time. Press ENTER to activate either the Date or the Time menu item, use the UP or DOWN arrow key to change the Date or the Time. After each change, the cursor moves to the right. When done, press ENTER to confirm and press MENU to exit. To change either the date or the time format, press the MENU key and select Time Settings > Date Fmt or Time Fmt. Press ENTER to activate either the Date Fmt or the Time Fmt menu item, use the UP or DOWN arrow key to change the Date or the Time format. When done, press ENTER to confirm and press MENU to exit.

The HW4 software can also be used to set the clock of the HygroLab C1 to the PC date and time.

Select which probe and/or parameters are shown on the display

Press the MENU key and select Display Settings > Mode. Press ENTER to activate the Mode menu item, use the UP or DOWN arrow key to select the display mode. Press ENTER to confirm and press MENU to exit.

Use the UP or DOWN arrow key to change the probe being displayed.

NOTE: The calculated parameter is shown only if enabled for the probe input that is selected (MENU > Probe Settings).

Practical advice for measuring humidity

The most common source of error when measuring relative humidity is a difference between the temperature of the probe and the temperature of the environment. At a humidity condition of 50 %RH, a temperature difference of 1°C (1.8 °F) typically results in an error of 3 %RH on relative humidity. It is always good practice to monitor the display for temperature stability. The probe should be given sufficient time to equilibrate with the environment to be measured. The larger the initial temperature difference between the probe and the environment to be measured, the more time temperature equilibration requires. This time can be shortened, and errors avoided, by using the probe configuration that fits best for your application. In extreme situations, condensation may occur on the sensors when the probe is colder than the environment. As long as the humidity / temperature limits of the humidity sensor are not exceeded, condensation does not alter the calibration of the sensor. However, the sensor has to dry out before it can provide a valid measurement. Non-moving air is an excellent insulator. When there is no air movement, surprising differences in temperature and humidity can be noted over short distances. Air movement at the probe generally results in measurements that are both faster and more accurate.

Aw Measurement

When enabled to measure water activity, the HygroLab C1 automatically displays humidity as Aw (1.000 Aw = 100 %RH) and offers the following options:

AwE mode:

The HygroLab C1 waits for the full equilibration of the measured product and probe. For most products, this takes from 30 to 60 minutes. The HygroLab C1 automatically detects equilibrium conditions (humidity and temperature) and ends the measurement at that time by freezing the display.

AwQuick mode:

The HygroLab C1 uses an algorithm to accelerate the water activity measurement and provides a result in typically 5 minutes. The measurement ends automatically and the display is frozen. When temperature conditions are stable (both at the product and probe), the value measured with the AwQuick mode is generally within ± 0.005 Aw of the value that would be obtained by waiting for full equilibration of the product and probe.

NOTE: See technical paper E-T-AW (Measuring Water Activity) for basic information on water activity and its applications.

Enabling Aw measurement

To enable the HygroLab C1 to measure water activity:

- Press the MENU key and select Aw Settings. Press ENTER to activate the menu.
- With the Enable Aw menu item highlighted, press ENTER and use the UP or DOWN arrow key to select ON. Press ENTER to confirm the selection.
- Optional: Use the DOWN arrow key to select the SaveResult menu item and press ENTER. Use the UP or DOWN arrow key to select either ON or OFF. Press ENTER to confirm the selection
- Use the DOWN arrow key to select the Mode menu item and press ENTER. Use the UP or DOWN arrow key to select either AwQuick or AwE. Press ENTER to confirm the selection.
- The settings for either the AwQuick or AwE function can be changed after using the UP or DOWN arrow key to highlight the setting and by pressing on ENTER. Use the UP or DOWN arrow key to change each digit. Press ENTER to move the cursor to the right. When done, press ENTER to save the value.
- Press MENU twice to fully exit the menu.

Using the AwQuick mode

The HygroLab C1 uses an algorithm to project the full equilibrium value (water activity) of the product sample:

- The value of the humidity signal is constantly monitored
- The stability of the temperature signal is constantly monitored
- After an initial period of time („dwell“ time), the humidity data is used to project the end value of the equilibration process (water activity). The measurement ends automatically as soon as the projected Aw value is stable. At that time, the HygroLab C1 freezes the display.


The measurement is automatically ended and typically requires about 5 to 6 minutes.

With the default „dwell“ time of 4 minutes, the difference between the AwQuick function and the conventional measurement method is typically 0.005 aw or less. The value of the „dwell“ time can be set by the user (see Aw Settings) and is a tradeoff between speed of measurement and accuracy. Generally, a longer „dwell“ time produces more accurate results but causes measurements to take longer.

The value of temperature shown at the end of the measurement is the average temperature during the measurement. The HygroLab C1 displays a trend indicator to the left of the temperature value. This is used to verify that temperature was stable during the measurement.

1. AWQ Reset:


The HygroLab C1 is ready to start measuring all connected probes

When ready to measure, press on the ENTER key.

2. AWQ Dwell:

The measurement is in the “dwell” phase.

AwQuick mode is in the dwell phase

Current humidity value


Current temperature value

Elapsed time, hh:mm:ss

Probe 1 is being displayed

3. AWQ Running:

At the end of the “dwell” phase the HygroLab C1 starts projecting the end result


Current Aw projection

Average temperature

4. AWQ Ended:

When the projection is stable for a probe being displayed, the HygroLab C1 automatically ends the measurement and freezes the display for that probe. A check mark appears to the left of the water activity value. In addition the HygroLab C1 can be set with the HW4 software to emit an acoustic signal (Beep) lasting 5 seconds when the last probe is done. The acoustic signal can be stopped by pushing a key (do not press the ENTER key). When several probes are connected to the HygroLab C1, use the UP or DOWN key to view each probe.


IMPORTANT: do not press the ENTER key until you have reviewed each probe.

5. AWQ Reset:

Write down the measurement for each probe or use the automatic data capture function. This function can be enabled from the keypad: MENU > Aw Settings > SaveResult > ON. The Aw value measured by probe 1 is saved to Data Bin 1. Probe 2 is saved to Data Bin 2, etc.

Press ENTER. The HygroLab C1 is ready to start a new measurement

Technical data (operation)

Number of probe inputs	4
Probe type:	HygroClip 2 probes
Operating limits at probe:	Depends on probe model
Operating limits at electronics:	-10...60 °C
	0...100 % RH, non condensing
Housing material:	ABS, Aluminum
Dimensions:	172 x 70 x 220mm
Weight:	1000 g

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ROTRONIC Instrument Corp, NY 11788, USA

Phone +1 631 427-3898, www.rotronic-usa.com

ROTRONIC South East Asia Pte Ltd, Singapore 339156

Phone +65 6294 6065, www.rotronic.com.sg

ROTRONIC Shanghai Rep. Office, Shanghai 200233, China

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ROTRONIC AG, CH-8303 Bassersdorf
ROTRONIC Messgeräte GmbH, D-76275 Ettlingen
ROTRONIC Instruments (UK) Ltd, Crompton Fields,